

Patient Name

Measured By

Date

Organisation

Purchase Order No.

Billing Address

Delivery Address

Custom Made

Pertex
Light

Pertex

GOLDPUNKT

1 Style and Fabric

Drug Tariff codes in red



Below Knee

Class 1
18-21mmHg☐
PX-21Class 2
23-32mmHg☐
PX-22Class 3
34-46mmHg☐
PX-23Class 2
23-32mmHg☐
GP-22Class 3
34-46mmHg☐
GP-23Class 4
49-70mmHg☐
GP-24

Thigh High

☐
PX-31☐
PX-32☐
PX-33☐
GP-32☐
GP-33☐
GP-34

MEASUREMENT TIPS

- Use a skin marker to record points
- Lengths taken as straight distances
- Measure up the limb starting with lengths
- Always take A-d measurement

= Circumference

= Straight lengths

2 Measurements

FO011
FO011-001 ver.10

g THIGH TOP
4cm below the inner gluteal fold

e AROUND KNEE

d BELOW KNEE
3cm below the popliteal crease

c WIDEST CALF

b ANKLE
3cm above the malleolus

h HEEL AROUND THE ANKLE
with the foot positioned in dorsi flexion

A FLOOR

a BASE OF TOES

HEEL TO BASE OF BIG TOE

HEEL TO TOP OF LONGEST TOE

Right Left

3 Toe



CLOSED

OPEN

☐☐

4 Colour



BEIGE

BLACK

☐☐

5 Fastening

Default position inside leg, between b and d



NONE

ZIP


VELCRO

☐☐☐

CM-92*

CM-93*

6 Grip Top



NONE

SILICONE

☐☐

CM-90*

7 Quantity

1 = 1 stocking

RIGHT

LEFT

COMMENTS e.g. alternate position for fastening

REPEAT SERIAL NO.

