



MEASUREMENTS

THIGH	g	Just below the groin
	f	Mid-thigh
	e ¹	Bottom of the thigh
	e ¹ -g	Follow the contours of the limb
KNEE	e ¹	Bottom of the thigh
	e	Over knee
	d	2 finger width below knee crease
	d-e ¹	Follow the contour of the limb
LEG	d	2 finger widths below knee crease
	c	Widest point of calf
	b ²	Mid-point between b ¹ and c
	b ¹	Point at which the achilles tendon changes into calf muscle
	b	Ankle
	b-d	Follow the contour of the limb
	h	Heel
FOOT	a ¹	Mid foot
	a	Base of Toes
	A-b	Straight distance
	a-heel	Straight distance

Patient Name	Order No.	Repeat Garment No.
Date Measured	Clinic / Hospital	
Measured by	Telephone	E-mail

1 QUANTITY

FOOT	Right	Left
Light EL-FC		
Strong ES-FC		



LEG	Right	Left
Light EL-LC		
Strong ES-LC		




KNEE ONLY	Right	Left
Light EL-KC		
Strong ES-KC		




THIGH ONLY	Right	Left
Light EL-TC		
Strong ES-TC		

2 LENGTH MEASUREMENTS Take measurements at back of leg



4 ADDITIONAL LINERS	STANDARD C circumference < 60cm		EXTRA WIDE C circumference > 60cm	
	PAIR OF A-d LINERS INCLUDED WHEN YOU ORDER FOOT & LEG	A-d	E-L-AD	E-L-ADW
	PAIR OF A-g LINERS INCLUDED WHEN YOU ORDER FOOT, LEG & THIGH	A-g	E-L-AG	E-L-AGW

3 CIRCUMFERENCE MEASUREMENTS



5 CHOOSE COLOUR	
SAND	BLACK

6 COMMENTS

CAUTION: Assess shape of limb and ensure sufficient padding is used to protect vulnerable areas. easywrap Strong fabric should not be used on ankles under 18cm. For ankle measurements under 18cm apply enough padding to bring the ankle measurement to above 18cm. For ankles between 18cm and 22cm easywrap Strong should only be applied after specialist referral and under supervision. For ankles under 18cm easywrap Light may be used after specialist referral and under supervision.