




Patient Name	PO No.	Repeat Garment No.
Date Measured	Clinic / DME	
Measured by	Telephone	E-mail

1

QUANTITY

FOOT	Right	Left
Light EL-FC		
Strong ES-FC		



LEG	Right	Left
Light EL-LC		
Strong ES-LC		



KNEE ONLY	Right	Left
Light EL-KC		
Strong ES-KC		



THIGH ONLY	Right	Left
Light EL-TC		
Strong ES-TC		

MEASUREMENTS GLOSSARY	
THIGH	g Just below the groin
	f Mid-thigh
	e¹ Bottom of the thigh
KNEE	e¹-g Follow the contours of the limb
	e¹ Bottom of the thigh
	e Over knee
	d 2 finger width below knee crease
	d-e¹ Follow the contour of the limb
LEG	d 2 finger widths below knee crease
	c Widest point of calf
	b² Mid-point between b ¹ and c
	b¹ Point at which the achilles tendon changes into calf muscle
	b Ankle
	b-d Follow the contour of the limb
	h Heel
FOOT	a¹ Mid foot
	a Base of Toes
	A-b Straight distance (floor to B landmark)
	a-heel Straight distance (base of toes to heel)

2

LENGTH MEASUREMENTS

Take measurements at back of leg

RIGHT	LEFT
g cm	g cm
e¹-g cm	e¹-g cm
e¹ cm	e¹ cm
d-e¹ cm	d-e¹ cm
d cm	d cm
b-d cm	b-d cm
A-b cm	A-b cm
b cm	b cm
a-heel cm	a-heel cm

3

CIRCUMFERENCE MEASUREMENTS

RIGHT	LEFT
g cm	g cm
f cm	f cm
e¹ cm	e¹ cm
e cm	e cm
d cm	d cm
c cm	c cm
b² cm	b² cm
b¹ cm	b¹ cm
b cm	b cm
h cm	h cm
a¹ cm	a¹ cm

4

ADDITIONAL LINERS

		STANDARD circumference < 60cm	EXTRA WIDE circumference > 60cm
PAIR OF A-d LINERS INCLUDED WHEN YOU ORDER LEG	A-d	E-L-AD	E-L-ADW
PAIR OF A-g LINERS INCLUDED WHEN YOU ORDER THIGH	A-g	E-L-AG	E-L-AGW

5

CHOOSE COLOR

SAND	BLACK
------	-------

6

COMMENTS

CAUTION: Assess shape of limb and ensure sufficient padding is used to protect vulnerable areas. easywrap Strong fabric should not be used on ankles under 18cm. For ankle measurements under 18cm apply enough padding to bring the ankle measurement to above 18cm. For ankles between 18cm and 22cm easywrap Strong should only be applied after specialist referral and under supervision. For ankles under 18cm easywrap Light may be used after specialist referral and under supervision.